

PRELOADED DMEK TISSUE

in your choice of three injectors

DMEK surgery is simple and fast with preloaded tissue prepared by Eversight, giving you more freedom in the OR.

Preloaded DMEK tissue service

Tissue for DMEK is prepared according to your specifications and available in one of three injectors. Using preloaded DMEK tissue streamlines DMEK surgery by saving you time and supplies in the operating room, allowing you to focus on patient care.

- DORC tube
- Modified Jones tube
- Micro Jones tube

Markings

- I-II, S or F marking

Graft sizes

- 7.0, 7.25, 7.5, 7.75, 8.0, 8.25 or 8.5mm

Tissue culturing


- A free-floating 5.0mm rim segment is included in the vial for culturing needs

DMEKConnect™


Simplify DMEK surgery with DMEKConnect.

DMEKConnect eliminates the need for an IV tubing attachment, reducing preparation time and eliminating steps in the OR.


Preloaded DMEK options



The **Modified Jones tube** is capped at the beveled end, DMEKConnect is attached at the back end.



A smaller incision option, the **Micro Jones tube** is also capped at the beveled end with a DMEKConnect attached at the back end.



The **DORC tube** is capped at the insertion end and connects easily to a syringe at the back end.

On-demand training

Eversight provides free online training to give you the resources you need, when you need them. To learn more, visit eversightvision.org/onlinetraining

*Preloaded DMEK is not available in the following countries: Austria, Belgium, Bulgaria, Croatia, Czech Republic, Denmark, Finland, France, Germany, Greece, Hungary, Ireland, Italy, Netherlands, Norway, Poland, Portugal, Romania, Serbia, Spain, Sweden, Switzerland, Turkey and United Kingdom.