

# TISSUE SERVICES

Eversight delivers exceptional customer service and access to the latest tissue processing techniques through a global network of community eye banks. All tissue provided by Eversight meets the highest quality standards established by the Eye Bank Association of America, U.S. Food and Drug Administration and Eversight Medical Advisory Committee.

## Tissue for transplant

Prepared according to surgeon preference.

### Cornea or whole eye

(penetrating or lamellar transplant; keratoprosthesis; keratolimbal allograft)

- Cornea with epithelium centration dot
- Cornea with trephine marking

### Cornea prepared for ALK or DSAEK

(Lamellar dissection by Moria microkeratome)

ALK

- DSAEK (90 $\mu$ -135 $\mu$ )
- Ultrathin DSAEK (40 $\mu$ -90 $\mu$ )

### Pre-loaded DSAEK

- 70-135 $\mu$

### Cornea prepared for DMEK

- Punched to size
- Pre-loaded

### Cornea prepared for LAK/IEK

Laser assisted keratoplasty  
Intralase enabled keratoplasty

### Cornea for glaucoma shunt coverage

- Stored in Optisol GS
  - » Whole
  - » Lamellar-prepared
  - » Bisected
- Long-term preserved in ethanol or glycerol
  - » Whole
  - » Bisected
  - » Anterior lenticule

### Sclera

- Whole
- Half
- Quarter
- Eighth
- 6mm disc

## Tissue for research

### Whole eye

### Posterior pole

### Cornea

### Conjunctiva

### Lens

### Limbal scleral rim remnant

### For a specific surgical service or tissue preparation not included here, contact:

Tissue Placement Department  
(866) 900-8119  
distribution@eversightvision.org

### For a specific research-related service, tissue preparation or preservation, contact:

Research Department  
(312) 469-5530  
research@eversightvision.org

*An Eversight representative is always available to accommodate emergency tissue requests.*