

SURGICAL PRACTICE SYSTEM

The Eversight Surgical Practice System features the Barron Artificial Anterior Chamber (AAC), providing a life-like, responsive alternative to whole globes for practicing surgical techniques.

Better corneal clarity vs. whole globes

• Outstanding visualization of anterior chamber depth

Versatile with exceptional control

- Multiple options for practicing surgical techniques including DSAEK, DMEK and preloaded.
- Easily adjust anterior and posterior pressure. Each system includes one full-thickness cornea.
- Upon request, the tissue can be delivered preassembled on the Barron AAC.

On-demand training

• Online ordering and instructional resources at eversightvision.org/sps.

Simplify DMEK surgery with **DMEKConnect**, the Eversight Modified Jones Tube and syringe connector. **DMEKConnect** eliminates the need for an IV tubing attachment, reducing preparation time and eliminating steps in the OR.

Standard Surgical Practice System

- Barron AAC
- Latex iris
- Seating block
- Tying forceps
- 15mL BSS
- 3mL syringe (2)
- 26g cannula
- Drape
- Gauze
- Surgical spears
- Surgical marker
- Shipping

International shipping requires additional fee.

DMEK Basic

- Standard Surgical Practice System, plus
- 3mL syringe
- Modified Jones Tube and **DMEKConnect**
- Donor cornea punch (7.5, 8 or 8.5mm)

DMEK Complete

- Standard Surgical Practice System, plus
- 3mL syringe
- Modified Jones Tube and **DMEKConnect**
- 30mL BSS
- Petri dish
- Surgical instruments: straight forceps, 15° microknife, 3mm cornea knife, crescent blade