

TISSUE SERVICES

Eversight delivers exceptional customer service and access to the latest tissue processing techniques through a global network of community eye banks. All tissue provided by Eversight meets the highest quality standards established by the Eye Bank Association of America, the U.S. Food and Drug Administration and the Eversight Medical Advisory Committee.

Tissue for transplant

Prepared according to surgeon preference.

Cornea or whole eye

(penetrating or lamellar transplant; keratoprosthesis; keratolimbal allograft)

- Cornea with epithelium centration dot

Cornea prepared for ALK or DSAEK*

(Lamellar dissection by Moria microkeratome)

- ALK
- DSAEK (90-135µm)
- Ultrathin DSAEK (40-90µm)

Preloaded DSAEK

- 70-135µm

*Other options available

Cornea prepared for DMEK

- Punched to size
- Preloaded Modified Jones Tube
- Preloaded DORC tube

Cornea prepared for LAK/IEK

- Laser-assisted keratoplasty
- Intralase-enabled keratoplasty

Cornea for glaucoma shunt coverage

- Stored in LIFE4°C
 - » Whole
 - » Lamellar-prepared
 - » Bisected
- Long-term preserved in ethanol or glycerol
 - » Whole
 - » Bisected
 - » Anterior lenticule

Sclera

- Whole

halo®

- Sterile pericardium
- Sterile sclera
- Sterile bisected/split thickness cornea
- Sterile full/split thickness cornea
- Sterile cornea (without sclera)
- Sterile cornea (with sclera)

Tissue for research

Service fees vary depending on researcher specified protocols.

Aqueous humor	Posterior pole
Conjunctiva	Retina
Iris	RPE/Choroid complex
Lens	Vitreous humor
Limbal scleral rim remnants	Whole eye
Optic nerve	

Contact

For a specific surgical service or tissue preparation not included here, contact:

Tissue Distribution
(866) 900-8119 | distribution@eversightvision.org

For a specific research-related service, tissue preparation or preservation, contact:

Research
(312) 469-5530 | research@eversightvision.org

An Eversight representative is always available to accommodate emergency tissue requests.